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(0.8 – 2.5 cm) thick.

Caution: Federal law (USA) restricts this device to sale by or on the order of a licensed practitioner.

Contraindications

- Do not use this device in an MRI environment, in an explosive atmosphere, or on infant or neonatal patients.
- This device is not defibrillation proof per IEC 60601-1:1990 clause 17h.

Warnings

- This device is intended only as an adjunct in patient assessment. It must be used in conjunction with other methods of assessing clinical signs and symptoms.
- The device must be able to measure the pulse properly to obtain an accurate SpO₂ measurement. Verify that nothing is hindering the pulse measurement before relying on the SpO₂ measurement.
- Operation of this device below the minimum amplitude of 0.3% modulation may cause inaccurate results.
- General operation of the device may be affected by the use of an electrosurgical unit (ESU).
- The use of accessories other than those specified in these instructions may result in increased electromagnetic emission and/or decreased immunity of this device.
- This device should not be used adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, the device should be observed carefully to verify normal operation.

Cautions

- This device is designed to determine the percentage of arterial oxygen saturation of functional hemoglobin. Significant levels of dysfunctional hemoglobin may affect measurement accuracy.
- This device has no audible alarms and is intended only for spot-checking.
- Inspect the sensor application site at least every 6 to 8 hours to ensure correct sensor alignment and skin integrity. Patient sensitivity to sensors may vary due to medical status or skin condition.
- Fluctuating or very bright light, moisture, blood pressure cuffs, infusion lines, venous pulsations, insufficient pulse signal, anemia, arterial catheters, nail polish, and/or artificial nails may degrade the device's performance.
- The device may not work when circulation is reduced. Warm or rub the finger, or re-position the device.
- This device's display will go blank after 30 seconds of no readings or poor readings.
- In some circumstances, the device may interpret motion as good pulse quality. Minimize patient motion as much as possible.
- Cardiac green and other intravascular dyes may affect the accuracy of the SpO₂ measurement.
- Do not autoclave or immerse this device in liquid or use caustic or abrasive cleaning agents.
- A flexible circuit connects the two halves. Do not twist or pull the flexible circuit or overextend the device's spring.
- Do not hang the lanyard from the device's flexible circuit/strain relief.
- A functional tester cannot be used to assess the accuracy of a pulse oximeter monitor or sensor.

Webmaster@the360.com Nonin produces the best pulse oximeters for the hospital and home care markets.

Nonin® 9600

People v. Murray

SA073164

Exh.

**MedAir™ LifeSense®
Capnograph and
Pulse Oximeter**

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